






VNG Medical Innovation System

VNG MEDICAL INNOVATION SYSTEM

 Plot No-197, JLPL Industrial Park
Sector-82, Mohali

 vngmedical@yahoo.in

 +91 87250 31075
0172 2970197

 www.vngmedical.com

Multipara Monitor Stand



Weight	1 kgs
Material	Aluminium
Height	300mm (Adj. Length)
Movement	Left/right, Up/Down, $\pm 10^\circ$
Basket	Optional
Basket material	SS/MS/Acrylic
Load capacity	20kgs



Basket material Stain Steel



Basket material Acrylic




Syringe Pump Stand

Weight	5 kgs
Height	1650mm
Base size	600mm
Material	SS/MS
Load capacity	25 kgs
Power cord	6 nos



Monitor Stand with Basket



Weight	5 kgs
Height	1250mm
Base size	600mm
Material	SS/MS
Base material	SS/Plastic
Basket optional	Single/Double
Basket material	MS Chrome plated/SS/Acrylic
Load capacity	25kgs

